



Patent Pending

COBBRA

Energized Cobb Elevator


- Effortless Dissection
- Minimizes Bleeding
- Safe Bone Preparation
- RF Tissue Preservation
- Reduces Operative Time
- Manual Cobb Ability



Patent Pending

COBBRA

Product Specifications and Technology Overview



Our Technology... See The Difference.

elliquance™ technology incorporates the patented Surgi-Max® Radiowave energy source which utilizes high frequency, low temperature radiowaves to perform traditional scalpel, scissor, electro-surgical and laser assisted procedures. The cell-specific tissue effect affords unparalleled surgical precision while sparing healthy tissue. The low temperature emission results in non-adherent bipolar performance which minimizes tissue trauma and eliminates frequent cleaning and instrument irrigation.



RF-COBB/M



RF-COBB/L

elliquance
Innovative Medical Solutions

2455 Grand Avenue, Baldwin, New York 11510 U.S.A.
Phone: (516) 277-9000 • Fax: (516) 277-9001 • www.elliquance.com SM 013 A