



# Surg- Pak

*4 Electrodes & 1 Handpiece*

- Convenient all-in-one kit
- Quick change electrodes
- Secure spin free electrodes
- Multiple procedure functionality
- 3 modality activation control
- Comfort rubber grip
- Ergonomic design



*Surg-e Pak is designed to be used exclusively with the elliquence family of generators*

# Product Specifications and Technology Overview



**SP-01**



## Our Technology... See The Difference.

elliquance™ technology incorporates the patented Surgi-Max® Radiowave energy source which utilizes high frequency, low temperature radiowaves to perform traditional scalpel, scissor, electrosurgical and laser assisted procedures. The cell-specific tissue effect affords unparalleled surgical precision while sparing healthy tissue. The low temperature emission results in non-adherent bipolar performance which minimizes tissue trauma and eliminates frequent cleaning and instrument irrigation.

	<p><b>ST-L5</b> 6mm Loop 4cm working length</p>
	<p><b>ST-B2</b> 2mm Ball 4cm working length</p>
	<p><b>ST-S</b> Spatula 4cm working length</p>
	<p><b>ST-W</b> Fine Wire 4cm working length</p>

**elliquance**  
Innovative Medical Solutions

2455 Grand Avenue, Baldwin, New York U.S.A.  
Phone: (516) 277-9000 • Fax: (516) 277-9001 • [www.elliquance.com](http://www.elliquance.com) SM 027 A