




RF Micro Fibre Kit

- Electrode length: 23 inches
- Ablation, Coagulation and Cutting Capabilites
- Monopolar Functionality
- Flexible
- 16 gauge insulated needle (included)



Product Specifications and Technology Overview





Our Technology... See The Difference.

elliquence™ technology incorporates the patented Surgi-Max® Radiowave energy source which utilizes high frequency, low temperature radiowaves to perform traditional scalpel, scissor, electrosurgical and laser assisted procedures. The cell-specific tissue effect affords unparalleled surgical precision while sparing healthy tissue. The low temperature emission results in non-adherent bipolar performance which minimizes tissue trauma and eliminates frequent cleaning and instrument irrigation.

