




Bi-Tip™


Micro Bipolar Handpiece

- Radiopaque
- Tactile Feedback
- Cut or Ablate
- 25cm Working Length
- Compatible with 4mm or larger working channels



Bi-Tip™

Product Specifications and Technology Overview



Our Technology... See The Difference.

elliquance™ technology incorporates the patented Surgi-Max® Radiowave energy source which utilizes high frequency, low temperature radiowaves to perform traditional scalpel, scissor, electro-surgical and laser assisted procedures. The cell-specific tissue effect affords unparalleled surgical precision while sparing healthy tissue. The low temperature emission results in non-adherent bipolar performance which minimizes tissue trauma and eliminates frequent cleaning and instrument irrigation.



BT-FT

elliquance
Innovative Medical Solutions

2455 Grand Avenue, Baldwin, New York 11510 U.S.A.
Phone: (516) 277-9000 • Fax: (516) 277-9001 • www.elliquance.com SM 007 A
Copyright 2008 © elliquance, LLC. All Rights Reserved.