



elliquence
Less Is More®

Trigger•Flex® Quad

- Compatible with Stenosis Spine Scopes
- Safe and Precise Nucleus Ablation and Spine Coagulation
- Helps provide a clear picture to determine anatomic landmarks
- Quick Control of bleeding to regain visualization
- Safe and Precise Annulus Modulation
- Steerable / Navigational Access
- Tactile Feedback / Probing
- Radiopaque



Technology Overview




Modes for Trigger-Flex® Quad with Surgi-Max®

Bipolar Hemo
Bipolar Turbo



Used exclusively with Surgi-Max®
Radiowave Technology




Articulated Surgical Precision
for targeted areas



Low Heat Temperatures
allow minimal damage
to surrounding tissue




Quad Bipolar Tip
for focused energy
in larger areas



Our Technology...See The elliquence Difference


elliquence has manufactured innovative radiofrequency systems and electrodes since 1969. The Trigger-Flex handpiece combined with the Surgi-Max High Frequency Radiowaves is designed for minimally invasive procedures such as endoscopic spine surgery. While laser and most electrosurgical devices are heat driven and imprecise, Surgi-Max provides unparalleled surgical precision at a much lower temperature allowing it to gently target tissue with minimal damage to surrounding healthy tissue. ^{A, B, C, D, E}



DTFQ

Working Length: 31cm

Outer Diameter: 5.5mm x 3.5mm



Warning and precautions:

- A) Hazardous electrical output. This equipment is for use only by qualified personnel.
- B) Only elliquence supplied and/or approved accessories should be used to ensure proper operation. The elliquence generator and its accessories are identified with the manufacturer information and reference identification number for proper traceability.
- C) Read the instruction manual in detail before use.
- D) Use of accessories or cables other than elliquence specified, as replacement parts for internal components, may result in increased emissions or decreased immunity of the elliquence RF generator.
- E) The coupling of any non-elliquence manufactured product, with an elliquence manufactured product invalidates all elliquence product warranties and all elliquence legal liabilities associated with such use.

elliquence
Less Is More®

2455 Grand Avenue • Baldwin, New York 11510-3531 U.S.A.
Phone: 516-277-9000 • Fax: 516-277-9001 • www.elliquence.com
Copyright 2016 © elliquence, LLC. All Rights Reserved
SM 049 A