



**elliquence**  
Less Is More®



**SURGI-MAX**  
**ULTRA**

## The Next Evolution in Medical Technology


- Unique Proprietary Ultra Wave Technology
- Enhanced cutting capabilities through different types of tissue (ie: muscle to fat)
- Minimal lateral heat spread that allows for shorter post op recovery
- Multiple modes provide versatile functionality
- 150W (Monopolar) / 170W (Bipolar) Output Power
- Cardio Floating (CF) approved

**Unique  
Ultra Wave  
Technology**

# SURGI-MAX<sup>®</sup> ULTRA

## Technology Overview

Product Code: IEC6-SU170



## Our Technology...See The Difference

In collaboration with doctors and surgeons from around the globe, elliquence has tailored it's focus to creating a unique energy source. The Surgi-Max Ultra has proprietary Ultra Wave technology that optimizes it's capabilities. This new technology has enhanced cutting capabilities through different types of tissue while minimizing the amount of lateral heat spread. The Surgi-Max Ultra offers multiple modalities to provide a versatile array of procedures and instrumentation.

### Monopolar Advantages

Proprietary Ultra Frequency can allow for simultaneous lower temperature tissue resection and coagulation.

- Wide array of electrode tips and sizes that allow for an wide variety of applications
- Multiple modes that allow for cutting, blended cutting and coagulation without carbonization
- Hemostasis with less heat
- Enhanced cutting capabilities through different types of tissue (ie: muscle to fat)
- Reduced lateral heat spread
- Optimal non-stick performance

### Bipolar Advantages

Proprietary Ultra Frequency can allow for bipolar coagulation with reduced carbonization.

- Optimal non-stick performance
- No irrigation required
- Reduction in lateral heat spread
- Reduced carbonization
- Improved hemostasis
- Improved safety
- Improved surgical workflow
- Less post op recovery

