

# MULT- PORTAL<sup>®</sup>



## What is Multi-e Portal?

A unique endoscopic approach that allows for two working channels to access the target pathology.

## MULT- PORTAL

Bridging the Gap Between  
Uniportal and Biportal Endoscopy

Uniportal Approach




Biportal Endoscopy

## What are the benefits of the Multi-e Portal approach when compared to Uniportal and Biportal Endoscopy?

- Multiple working ports enhance flexibility and accessibility
- Minimizes tissue and multifidus muscle dissection
- Enables creation of a surgical cavity from the primary port
- Allows docking a cannula on bone to reduce fatigue
- Maintains irrigation flow with an independent irrigation port
- Reduces irrigation volume and avoids relying on the working port for fluid management
- Ideal for wide or multi-level decompression
- Potential to reduce surgical time

2 Working Channel Ports



Scan Here to Schedule Your Training

 elliquence  
Less Is More<sup>®</sup>

	QTY	DESCRIPTION
	1	15-2906 - Short Working Port, Half Port, OD: 12.0mm, ID: 10.0mm, WL:60.0mm
	1	11-2308 - Coaxial Step Dilator (IL/ST #1) - OD: 3.95mm, WL: 240mm 11-2309 - Coaxial Step Dilator (IL/ST #2) - OD: 6.8mm, ID:4.0mm, WL: 220mm
	1	15-2440 - Chisel, 0°, OD: 5.0mm, WL: 148.2mm
	1	15-2534 - Curette, 35°, OD: 3.0mm, WL: 150.0mm
	1	15-2535 - Curette, 35°, OD: 5.0mm, WL: 150.0mm
	1	15-2401 - Freer Elevator, Double Ended 15°/30°, OD: 4.0mm, WL: 118.5mm
	1	15-2404 - Nerve Probe, L-Shape, Shaft OD: 4.0mm, Tip OD: 1.5mm, WL: 152mm
	1	15-2310 - Muscle Dissector, OD: 12.0mm, WL: 122.0mm
	1	11-2511 - Mallet
	1	15-1641 - Cup Forceps, OD: 4.0mm, WL: 158mm
	1	15-1642 - Cup Forceps, Upward Angle, OD: 4.0mm, WL:158mm
	1	15-1934 - Kerrison, 90°, OD: 3.0mm, WL: 150.0mm
	1	15-1939 - Kerrison, Bent Upwards, 130°, OD: 3.0mm, WL: 150.0mm
	2	12-1947 - Kerrison Handle (Includes 2)

**MULT- PORTAL** Compatible Instruments Not Included in the Mult-e Portal Set

	DESCRIPTION
	15-2321 - Retractor, Straight, Tapered, 62°
	15-2442 - Chisel, 20°, OD: 5.0mm, WL: 148.2mm 15-2444 - Chisel, 38°, OD: 5.0mm, WL: 148.2mm 15-2449 - Chisel, Hockey 35°, OD: 5.0mm, WL: 148.2mm
	15-2532 - Curette, 45°, OD: 3.0mm, WL: 150.0mm 15-2533 - Curette, 130°, OD: 3.0mm, WL: 120.0mm 15-2536 - Curette, 90°, OD: 5.0mm, WL: 150.0mm
	15-1924 - Kerrison, 90°, OD: 2.0mm, WL: 150mm
	15-2920 - Suction Tube, 145°, OD: 3.5mm, ID: 2.4mm, WL: 145.4mm



Scan Here to Schedule Your Training